

FAST SOLUTIONS FOR EFFECTIVE DISINFECTION WITHIN 1 MINUTE

SPORICIDAL • TUBERCULOCIDAL • BACTERICIDAL • VIRUCIDAL • FUNGICIDAL

- Unique patented formulation eliminates spores, viruses, bacteria and fungi (including C. Difficile, Hepatitis B/C, HIV, H1N1, MRSA, TB and NDM-1) **within 1 minute**
- Non-abrasive, non-corrosive and non-volatile making it safe to use on all surfaces
- Pleasant scent

| Activity Claim | Contact Time | Test Protocol | Test Organism |
|-------------------------------------|--------------|--------------------|---|
| Bactericidal | 1 Minute | EN1276 | Escherichia coli, Staphylococcus aureus, Pseudomonas aeruginosa, Salmonella typhirium, Enterobacter cloacae, Enterococcus hirae, Lactobacillus plantarum, Klebsiella pneumonia, Enterococcus faecalis, Proteus mirabilis, Streptococcus pyogenes, Salmonella choleraesuis, Listeria monocytogenes |
| | | EN13727 | MRSA, Klebsiella pneumonia, Staphylococcus aureus, Pseudomonas aeruginosa, Enterococcus hirae, Acinetobacter Baumannii |
| | | EN13697 | Pseudomonas aeruginosa, Staphylococcus aureus, Escherichia coli, Enterococcus hirae |
| Fungicidal | 1 Minute | EN1650 | Candida albicans, Aspergillus niger, Tricophyton mentagrophytes |
| | | EN13624 | Aspergillus niger, Candida albicans |
| | | EN1675 | Aspergillus niger, Candida albicans |
| Tuberculocidal/ Mycobactericidal | 1 Minute | EN14348 EN14204 | Mycobacterium terrae, Mycobacterium avium Mycobacterium avium |
| Virucidal | 1 Minute | EN14476 | HIV, Human Corona virus, Human Influenza virus, Herpes simplex virus, Swine Influenza A virus |
| | | DW/RKI | Vaccinia virus, BVDV (Hepatis C surrogate) |
| Sporicidal | 1 Minute | EN13704 | Clostridium difficile spores, Bacillus subtilis spores, Clostridium sporogenes spores, bacillus cereus spores |

Meets EU requirements including: EN1276; EN1650; EN14348; EN13704; EN14476; EN13727; EN1656; EN13624; EN1657; EN14204

VAH/DGHM listed as a surface disinfectant, tested according to new standard method.